

# Flexmobil i6 Comfort II

With the Flexmobil i6 Comfort II, we increase the sitting comfort one step more, for those who sit frequently and for a long time in their electric wheelchairs. The Flexmobil i6 Comfort II is a powered wheelchair for indoor use that provides a new freedom of movement at home and at work. The chair is equipped with an innovative chassis with a more generous seat height and lower noise level. The new Comfort II represents an evolution of the previous Comfort chair, where a new sliding back support function creates a higher level of sitting comfort through the introduction of new cushions. In addition, it's also easier for the user to adjust the chair.

Type of impairment

minor  moderate  extensive

indoor  outdoor  child  adult



## SEATING SYSTEMS

### The Comfort II



The comfortable Comfort seat provides ergonomic support for all of the body, allowing you to sit well for the whole day. The seat is available in three widths and with a low or high back. You can also choose either Base, which is built from several foam layers, or Contour, which combines multiple foam layers with other elements made from memory foam. The seat position is fully adjustable up to 60 cm, allowing long-legged users to find a comfortable sitting position also.

Seat (width x depth)	40x40, 40x45, 40x50, 40x55, 45x40, 45x45, 45x50, 45x55, 50x40, 50x45, 50x50, 50x55 cm
Seat width (between armrests)	35-70 cm
Seat depth	40-60 cm
Seat height (height to underneath the seat)	39-64 cm, 46-86 cm
Seat tilt, electrical, crank & turnbuckle	-19° to 20°
Back (width x height)	40x45, 45x45, 50x45, 40x54, 50x54, 50x54 cm
Adjustable back angle	20°
Adjustable back support	47-60 cm
Upholstery fabric	Contour: Antracit Trikåby / Black Baltex Base: Black Baltex
Armrest, height, adjustable	10-39 cm
Armrest, platform	30x8x3 cm

The seat height is adjusted electrically and is available in high and low versions. Both are listed after each other in the technical data in the order "high, low".

## CHASSIS

### The Flexmobil i6



Flexmobil i6 makes it easy to both "sit near the floor" or "stand up and work" thanks to the generous min and max positions. The chassis is available with two different actuators. Thanks to the central placement of the drive wheel, the Flexmobil i6 is very easy to drive. The chair turns on its own axis, with an extremely small turning circle that doesn't need much space.

Chassis, width	57 cm
Chassis, length	80 cm
Total weight	110 kg
User weight, max.	150 kg
Turning circle	96 cm
Max speed	4,5 km/h
Max driving time	8,5 hr
Electronics	Dynamic/PG
Obstruction clearance	3 cm

Accessory code E09