

Junior

The Junior chair has been developed to help children/juniors with reduced mobility maintain a good posture during all the different activities of the day at school or nursery, at meal times and while playing. A broad range of accessories and easily managed functional settings make it easy to adapt the chair to the needs of the individual child. Junior has its own Junior52 chassis, but is also available with a Flexbase chassis.

Type of impairment

minor moderate extensive

indoor outdoor child adult



SEATING SYSTEMS

SitRite



The SitRite's advantages are hidden under the surface of the seat, in the two foam materials with different densities. The lower, harder, foam section helps to hold the pelvis in a good position, while the upper, softer, foam section achieves a distribution of pressure by being thicker under the seat bones than it is under the thighs. The SitRite is an excellent example of how design and density interact to create better balance, straightening of the spine and control of hand movements. All of which combine to provide the best levels of comfort and stability.

The SitRite's back has an ergonomic form, which supports a large part of the trunk without impairing the arm's freedom of movement. Pressure is distributed over a large area, and the trunk is also supported at the sides, which enhances comfort, increases stability and reduces tiredness when the user has reduced sitting capability. When the trunk is stable, the control of the hands and head's movements is also improved.

UPHOLSTERY FABRIC



Black Skai



Black fabric

Seat (width x depth)	29x32, 32x36, 36x40, 36x45, 40x40, 40x45 cm
Seat width (between armrests)	27-52 cm
Seat depth	32-45 cm
Seat height Gas (height to underneath the seat)	41-56 cm, 48-68 cm
Seat height Electrical (height to underneath the seat)	41-58 cm, 45-70 cm
Seat tilt electrical & gas	-20° to 12°
Seat tilt crank & turnbuckle	-16° to 20°
Back support, height, adjustable	37-54 cm
Back (width x height)	30x37, 37x43 cm
Adjustable back angle	30°
Upholstery fabric	Black Atlantic
Armrest, height, adjustable	0-26 cm
Armrest, platform	25x8x3 cm

The seat height is adjusted with the help of gas or electricity, where both versions are available in a high or low seat height version. Both are listed after each other in the technical data in the order "high, low".

CHASSIS

Junior52



The Junior52 chassis is very similar to its Flexbase equivalent, but it is narrower and more compact in order to enable it to fit in environments where space is at a premium, such as classrooms. However, Junior52 still offers plenty of kick space, and is very easy to roll even on surfaces that are not so smooth. The brake is easily accessible without impeding lateral movement. A rear wheel brake is fitted as standard, but it is easy to change to front wheel or four-wheel brakes. The two piece chassis means the chair remains level even on surfaces that are not so smooth.

Chassis, width	52 cm
Chassis, length	55 cm
Total weight	22 kg (gas), 29 kg (electrical)
User weight, max.	80 kg
Brake	Chassis mounted
Wheel size	100 mm

Accessory code E04