

Miniflex SitRite

The Miniflex SitRite is an indoor powered wheelchair with minimal external dimensions and a simple reliable design. Its small versatile format is extremely user friendly as is the advanced work chair seat, which has been developed jointly with experts in seat ergonomics. It allows the user to access tables and work surfaces with ease. Miniflex drives the front wheels and the slightly larger drive wheels allow you to drive over thresholds with ease. The jointed rear wheel pair means that the wheelchair always remains in contact with the floor.

Type of impairment

minor moderate extensive

indoor outdoor child adult



SEATING SYSTEMS

The SitRite



The SitRite's advantages are hidden under the surface of the seat, in the two foam materials with different densities. The lower, harder, foam section helps to hold the pelvis in a good position, while the upper, softer, foam section achieves a distribution of pressure by being thicker under the seat bones than it is under the thighs. The SitRite is an excellent example of how design and density interact to create better balance, straightening of the spine and control of hand movements. All of which combine to provide the best levels of comfort and stability.

The SitRite's back has an ergonomic form, which supports a large part of the trunk without impairing the arm's freedom of movement. Pressure is distributed over a large area, and the trunk is also supported at the sides, which enhances comfort, increases stability and reduces tiredness when the user has reduced sitting capability. When the trunk is stable, the control of the hands and head's movements is also improved.

UPHOLSTERY FABRIC



Black
Skai



Black
fabric

Seat (width x depth)	40x40, 40x45, 45x45, 45x50, 50x50 cm
Seat width (between armrests)	38-56 cm
Seat depth	40-59 cm
Seat height (height to underneath the seat)	40-58, 45-70 cm
Seat tilt, electrical, crank & turnbuckle	-22° to 15°
Back support, height, adjustable	37-62 cm
Back (width x height)	37x43, 42x45, 47x47 cm
Adjustable back angle	30°
Upholstery fabric	Black Atlantic
Armrest, height, adjustable	0-26 cm
Armrest, platform	30x8x3 cm

The seat height is adjusted electrically and is available in high and low versions. Both are listed after each other in the technical data in the order "high, low".

CHASSIS

The Miniflex



The market's most compact chair chassis for indoor use is no larger than a traditional mobility chair. The easily manoeuvred and compact design makes it easy to get close to tables or work surfaces. The chassis has front-wheel drive with a jointed rear axle, which ensures that the rear wheel pair is always in contact with the floor.

Chassis, width	57 cm
Chassis, length	72 cm
Total weight	60 kg
User weight, max.	130 kg
Turning circle	120 cm
Max speed	4 km/h
Max distance	approx. 15 km
Electronics	Dynamic/PG
Obstruction clearance	3 cm
Battery	12V/18Ah

Accessory code E10