

Mini

Mini is a safe, well-balanced small rollator that is equally at home both indoors and outdoors. It has been designed with children in mind, but can also be used by adults who are on the shorter side. At only 53 cm wide, the Mini is easy to operate.

The seat is stable. Comfort Frame™ technology allows the rollator to glide smoothly over uneven surfaces inside and out.

Its cross-function allows Mini to be easily folded up and stored standing.

The convenient handle height adjustment feature means that Mini can be tailored to each individual user. The wires are built into the handle, preventing them from catching on surrounding objects. Mini has a patented braking function (VABS), consisting of a "smart" linkage system. Both the service and parking brakes are ergonomically designed and easily accessible, giving the rollator a high level of safety.

A wide variety of accessories is available for the Mini, enabling it to satisfy the needs of a diverse range of users. One of these accessories is an articulation brake.

MODEL STANDARD

Mini

Weight	7.5 kg
User weight	125 kg
Width	53 cm
Length	74 cm
Seat width	36 cm
Height	68-84 cm
Seat height	47-51 cm
Width when folded	27 cm
Wheel dimensions	20x3.2 cm

MODEL LOW

Mini

Weight	7.1 kg
User weight	125 kg
Width	53 cm
Length	74 cm
Seat width	36 cm
Height	63-75 cm
Seat height	46 cm
Width when folded	27 cm
Wheel dimensions	20x3.2 cm

- indoor
- outdoor
- child
- adult

Children's
rollator for
both indoor and
outdoor use

