

Volaris HP

Volaris HP provides mobility support and is designed to be a walking aid for people suffering from weakness on one side of the body (Hemiplegia). The product can be used both indoors and outdoors.

The handle can be adjusted in multiple dimensions; horizontally, vertically, depth, and angle. For users who want to have their arms in a cupped position, the handle can also be equipped with forearm supports, all with the same multi-adjustment capability. The extra soft and wide wheels absorb and reduce vibrations on every step. The height of the handles and seat can be easily adjusted, meaning that the rollator is adaptable to all users' needs.

The signature Volaris cross design allows the rollator to travel smoothly across uneven surfaces and stand stably on all four wheels, including when it is folded up for storage or during transport.

Volaris HP features an extra wide, 48 cm seat that provides a stable surface to sit on. The maximum user weight of the rollator is 160 kg.

STANDARD MODEL HP

Weight	9.6 kg
User weight (maximum)	160 kg
Width	66 cm
Length	74 cm
Seat width	48 cm
Handle height, brake control	76-102 cm
Handle height, HP handle	92-118 cm
Handle height with forearm support	96-122 cm
Seat height	48.5-56 cm
Width when folded	30 cm
Wheel dimensions	20x4.2 cm

ITEM NUMBER

Volaris HP Left	1460100HPL
Volaris HP Right	1460100HPR
Forearm support for HP (accessory)	2466115

☒ indoor ☒ outdoor ☐ child ☒ adult

